

GSDMB Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51638-1**

Capture Antibody Information

Catalog Number: 67527-2-PBS	Clone ID: 1B12G9	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: gasdermin B
Isotype: IgG1	GenBank: BC025682	Gene ID: 55876
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag28332	

Detection Antibody Information

Catalog Number: 67527-3-PBS	Clone ID: 2G2E7	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: gasdermin B
Isotype: IgG1	GenBank: BC025682	Gene ID: 55876
Purification Method: Protein G purification	Immunogen Catalog Number: Ag28332	

Applications

Tested Applications: Cytometric bead array	Range: 3.125-400 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	--

Product Information

MP51638-1 targets GSDMB in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: GSDMB Monoclonal antibody, PBS Only (Capture) 67527-2-PBS (1B12G9). 100 μ g. Concentration 1 mg/mL.

Detection antibody: GSDMB Monoclonal antibody, PBS Only (Detector) 67527-3-PBS (2G2E7). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

