

For Research Use Only

# TAF6L Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51658-1

## Capture Antibody Information

Catalog Number:  
67569-2-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G Magarose purification

Clone ID:  
2F3H9

Reactivity:  
human

GenBank:  
BC008785

Immunogen Catalog Number:  
Ag8202

Conjugate:  
Unconjugated

Full name:  
TAF6-like RNA polymerase II, p300/CBP-associated factor (PCAF)-associated factor, 65kDa

Gene ID:  
10629

## Detection Antibody Information

Catalog Number:  
67569-3-PBS

Host:  
Mouse

Isotype:  
IgG2a

Purification Method:  
Protein A Magarose purification

Clone ID:  
1F10E9

Reactivity:  
human

GenBank:  
BC008785

Immunogen Catalog Number:  
Ag8202

Conjugate:  
Unconjugated

Full name:  
TAF6-like RNA polymerase II, p300/CBP-associated factor (PCAF)-associated factor, 65kDa

Gene ID:  
10629

## Applications

Tested Applications:  
Cytometric bead array

Range:  
1.563-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP51658-1 targets TAF6L in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TAF6L Monoclonal antibody, PBS Only (Capture) 67569-2-PBS (2F3H9). 100  $\mu$ g. Concentration 1 mg/mL.

Detection antibody: TAF6L Monoclonal antibody, PBS Only (Detector) 67569-3-PBS (1F10E9). 100  $\mu$ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array sample test of MP51658-1, TAF6L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67569-2-PBS. Detection antibody: 67569-3-PBS.



Cytometric bead array standard curve of MP51658-1, TAF6L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67569-2-PBS. Detection antibody: 67569-3-PBS. Standard: Ag8180. Range: 1.563-200 ng/mL.