

USP15 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP51664-1

Capture Antibody Information

Catalog Number:
67557-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
1C9G3
Reactivity:
human, mouse, rat
GenBank:
BC020688
Immunogen Catalog Number:
Ag5712

Conjugate:
Unconjugated
Full name:
ubiquitin specific peptidase 15
Gene ID:
9958

Detection Antibody Information

Catalog Number:
67557-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1C3C2
Reactivity:
human
GenBank:
BC020688
Immunogen Catalog Number:
Ag5712

Conjugate:
Unconjugated
Full name:
ubiquitin specific peptidase 15
Gene ID:
9958

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.


Product Information

MP51664-1 targets USP15 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: USP15 Monoclonal antibody, PBS Only 67557-1-PBS (1C9G3). 100 µg. Concentration 1 mg/mL.
Detection antibody: USP15 Monoclonal antibody, PBS Only (Detector) 67557-2-PBS (1C3C2). 100 µg. Concentration 1 mg/mL.
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.
Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51664-1, USP15 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67557-1-PBS. Detection antibody: 67557-2-PBS.



Cytometric bead array standard curve of MP51664-1, USP15 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67557-1-PBS. Detection antibody: 67557-2-PBS. Standard: Ag5712. Range: 0.781-100 ng/mL.