

HMGN2 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51668-1**

Capture Antibody Information

Catalog Number: 61074-1-PBS	Clone ID: 1B9D6	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: high-mobility group nucleosomal binding domain 2
Isotype: IgG1	GenBank: BC014644	Gene ID: 3151
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag28531	

Detection Antibody Information

Catalog Number: 61074-2-PBS	Clone ID: 2D5D3	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: high-mobility group nucleosomal binding domain 2
Isotype: IgG1	GenBank: BC014644	Gene ID: 3151
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag28531	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP51668-1 targets HMGN2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: HMGN2 Monoclonal antibody, PBS Only (Capture) 61074-1-PBS (1B9D6). 100 µg. Concentration 1 mg/mL.

Detection antibody: HMGN2 Monoclonal antibody, PBS Only (Detector) 61074-2-PBS (2D5D3). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51668-1, HMG2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61074-1-PBS. Detection antibody: 61074-2-PBS.



Cytometric bead array standard curve of MP51668-1, HMG2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61074-1-PBS. Detection antibody: 61074-2-PBS. Standard:Ag28531. Range: 0.781-100 ng/mL.