

For Research Use Only

VANGL2 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51674-1

Capture Antibody Information

Catalog Number:
67273-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
2B5E3

Reactivity:
human

GenBank:
BC103920

Immunogen Catalog Number:
Ag16980

Conjugate:
Unconjugated

Full name:
vang-like 2 (van gogh, Drosophila)

Gene ID:
57216

Detection Antibody Information

Catalog Number:
67273-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
2H8G1

Reactivity:
human

GenBank:
BC103920

Immunogen Catalog Number:
Ag16980

Conjugate:
Unconjugated

Full name:
vang-like 2 (van gogh, Drosophila)

Gene ID:
57216

Applications

Tested Applications:
Cytometric bead array

Range:
3.125-400 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51674-1 targets VANGL2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: VANGL2 Monoclonal antibody, PBS Only (Capture) 67273-2-PBS (2B5E3). 100 μ g. Concentration 1 mg/mL.

Detection antibody: VANGL2 Monoclonal antibody, PBS Only (Detector) 67273-3-PBS (2H8G1). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51674-1, VANGL2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67273-2-PBS. Detection antibody: 67273-3-PBS.



Cytometric bead array standard curve of MP51674-1, VANGL2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67273-2-PBS. Detection antibody: 67273-3-PBS. Standard: Ag16980. Range: 3.125-400 ng/mL.