

ELAVL4 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: MP51680-1

Capture Antibody Information

Catalog Number: 67835-2-PBS	Clone ID: 3E8C12	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: ELAV (embryonic lethal, abnormal vision, Drosophila)-like 4 (Hu antigen D)
Isotype: IgG1	GenBank: BC036071	Gene ID: 1996
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag3678	

Detection Antibody Information

Catalog Number: 67835-3-PBS	Clone ID: 2B1B6	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: ELAV (embryonic lethal, abnormal vision, Drosophila)-like 4 (Hu antigen D)
Isotype: IgG2a	GenBank: BC036071	Gene ID: 1996
Purification Method: Protein A Magarose purification	Immunogen Catalog Number: Ag3678	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP51680-1 targets ELAVL4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ELAVL4 Monoclonal antibody, PBS Only (Capture) 67835-2-PBS (3E8C12). 100 µg. Concentration 1 mg/mL.

Detection antibody: ELAVL4 Monoclonal antibody, PBS Only (Detector) 67835-3-PBS (2B1B6). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51680-1, ELAVL4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67835-2-PBS. Detection antibody: 67835-3-PBS.



Cytometric bead array standard curve of MP51680-1, ELAVL4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67835-2-PBS. Detection antibody: 67835-3-PBS. Standard:Ag3678. Range: 0.781-100 ng/mL