

FBXW7 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51687-1**

Capture Antibody Information

Catalog Number: 67234-2-PBS	Clone ID: 1F10C6	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: F-box and WD repeat domain containing 7
Isotype: IgG1	GenBank: NM_001349798	Gene ID: 55294
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag28471	

Detection Antibody Information

Catalog Number: 67234-1-PBS	Clone ID: 2A10E5	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: F-box and WD repeat domain containing 7
Isotype: IgG1	GenBank: NM_001349798	Gene ID: 55294
Purification Method: Protein A purification	Immunogen Catalog Number: Ag28471	

Applications

Tested Applications: Cytometric bead array	Range: 0.391-50 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP51687-1 targets FBXW7 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FBXW7 Monoclonal antibody, PBS Only (Capture) 67234-2-PBS (1F10C6). 100 µg. Concentration 1 mg/mL.

Detection antibody: FBXW7 Monoclonal antibody, PBS Only (Detector) 67234-1-PBS (2A10E5). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51687-1, FBXW7 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67234-2-PBS. Detection antibody: 67234-1-PBS.



Cytometric bead array standard curve of MP51687-1, FBXW7 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67234-2-PBS. Detection antibody: 67234-1-PBS. Standard: Ag28471. Range: 0.391-50 ng/mL.