

SNX5 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51693-1**

Capture Antibody Information

Catalog Number:
67665-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1A2A1

Reactivity:
human

GenBank:
BC093623

Immunogen Catalog Number:
Ag12330

Conjugate:
Unconjugated

Full name:
sorting nexin 5

Gene ID:
27131

Detection Antibody Information

Catalog Number:
67665-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
2A1C9

Reactivity:
human

GenBank:
BC093623

Immunogen Catalog Number:
Ag12330

Conjugate:
Unconjugated

Full name:
sorting nexin 5

Gene ID:
27131

Applications

Tested Applications:
Cytometric bead array

Range:
0.391-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51693-1 targets SNX5 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SNX5 Monoclonal antibody, PBS Only (Capture) 67665-2-PBS (1A2A1). 100 µg. Concentration 1 mg/mL.

Detection antibody: SNX5 Monoclonal antibody, PBS Only (Detector) 67665-3-PBS (2A1C9). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51693-1, SNX5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67665-2-PBS. Detection antibody: 67665-3-PBS.



Cytometric bead array standard curve of MP51693-1, SNX5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67665-2-PBS. Detection antibody: 67665-3-PBS. Standard:Ag12330. Range: 0.391-50 ng/mL.