

For Research Use Only

POLE Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51704-1

Capture Antibody Information

Catalog Number:
67526-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
3G4A3

Reactivity:
human

GenBank:
BC021559

Immunogen Catalog Number:
Ag30055

Conjugate:
Unconjugated

Full name:
polymerase (DNA directed), epsilon

Gene ID:
5426

Detection Antibody Information

Catalog Number:
67526-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
2A12G4

Reactivity:
human

GenBank:
BC021559

Immunogen Catalog Number:
Ag30055

Conjugate:
Unconjugated

Full name:
polymerase (DNA directed), epsilon

Gene ID:
5426

Applications

Tested Applications:
Cytometric bead array

Range:
3.125-400 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51704-1 targets POLE in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: POLE Monoclonal antibody, PBS Only (Capture) 67526-2-PBS (3G4A3). 100 μ g. Concentration 1 mg/mL.

Detection antibody: POLE Monoclonal antibody, PBS Only (Detector) 67526-3-PBS (2A12G4). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

