

For Research Use Only

SLC36A3 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51713-1

Capture Antibody Information

Catalog Number:
67929-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2C1D1
Reactivity:
human
GenBank:
BC101092
Immunogen Catalog Number:
Ag16556

Conjugate:
Unconjugated
Full name:
solute carrier family 36
(proton/amino acid symporter),
member 3
Gene ID:
285641

Detection Antibody Information

Catalog Number:
67929-3-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
1H9D7
Reactivity:
human
GenBank:
BC101092
Immunogen Catalog Number:
Ag16556

Conjugate:
Unconjugated
Full name:
solute carrier family 36
(proton/amino acid symporter),
member 3
Gene ID:
285641

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51713-1 targets SLC36A3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: SLC36A3 Monoclonal antibody, PBS Only (Capture) 67929-2-PBS (2C1D1). 100 μ g. Concentration 1 mg/mL.
Detection antibody: SLC36A3 Monoclonal antibody, PBS Only (Detector) 67929-3-PBS (1H9D7). 100 μ g. Concentration 1 mg/mL.
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.
Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

