

HIP1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51715-1**

Capture Antibody Information

Catalog Number: 67163-2-PBS	Clone ID: 1F8C2	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: huntingtin interacting protein 1
Isotype: IgG1	GenBank: BC110545	Gene ID: 3092
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag17982	

Detection Antibody Information

Catalog Number: 67163-3-PBS	Clone ID: 1H1C1	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: huntingtin interacting protein 1
Isotype: IgG1	GenBank: BC110545	Gene ID: 3092
Purification Method: Protein G purification	Immunogen Catalog Number: Ag17982	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	--

Product Information

MP51715-1 targets HIP1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: HIP1 Monoclonal antibody, PBS Only (Capture) 67163-2-PBS (1F8C2). 100 µg. Concentration 1 mg/mL.

Detection antibody: HIP1 Monoclonal antibody, PBS Only (Detector) 67163-3-PBS (1H1C1). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51715-1, HIP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67163-2-PBS. Detection antibody: 67163-3-PBS.



Cytometric bead array standard curve of MP51715-1, HIP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67163-2-PBS. Detection antibody: 67163-3-PBS. Standard: Ag17982. Range: 0.781-100 ng/mL.