

For Research Use Only

RAB11B Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51735-1

Capture Antibody Information

Catalog Number:
67780-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
1G7D6
Reactivity:
human, mouse, rat
GenBank:
BC110081
Immunogen Catalog Number:
Ag25429

Conjugate:
Unconjugated
Full name:
RAB11B, member RAS oncogene family
Gene ID:
9230

Detection Antibody Information

Catalog Number:
67780-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
2D2B1
Reactivity:
human
GenBank:
BC110081
Immunogen Catalog Number:
Ag25429

Conjugate:
Unconjugated
Full name:
RAB11B, member RAS oncogene family
Gene ID:
9230

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51735-1 targets RAB11B in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RAB11B Monoclonal antibody, PBS Only (Capture) 67780-1-PBS (1G7D6). 100 µg. Concentration 1 mg/mL.

Detection antibody: RAB11B Monoclonal antibody, PBS Only (Detector) 67780-2-PBS (2D2B1). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51735-1, RAB11B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67780-1-PBS. Detection antibody: 67780-2-PBS.



Cytometric bead array standard curve of MP51735-1, RAB11B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67780-1-PBS. Detection antibody: 67780-2-PBS. Standard: Ag25429. Range: 1.563-200 ng/mL