

CBL Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51740-1**

Capture Antibody Information

Catalog Number: 66576-2-PBS	Clone ID: 2G11H5	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: Cas-Br-M (murine) ecotropic retroviral transforming sequence
Isotype: IgG1	GenBank: BC132733	Gene ID: 867
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag23033	

Detection Antibody Information

Catalog Number: 66576-3-PBS	Clone ID: 2H5C4	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: Cas-Br-M (murine) ecotropic retroviral transforming sequence
Isotype: IgG1	GenBank: BC132733	Gene ID: 867
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag23033	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP51740-1 targets CBL in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CBL Monoclonal antibody, PBS Only (Capture) 66576-2-PBS (2G11H5). 100 µg. Concentration 1 mg/mL.

Detection antibody: CBL Monoclonal antibody, PBS Only (Detector) 66576-3-PBS (2H5C4). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

