

MYL6 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51759-1**

Capture Antibody Information

Catalog Number:
68142-1-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
2G4C9

Reactivity:
human, mouse, rat, pig, rabbit

GenBank:
BC006781

Immunogen Catalog Number:
Ag13091

Conjugate:
Unconjugated

Full name:
myosin, light chain 6, alkali, smooth muscle and non-muscle

Gene ID:
4637

Detection Antibody Information

Catalog Number:
68142-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
2B3D6

Reactivity:
human

GenBank:
BC006781

Immunogen Catalog Number:
Ag13091

Conjugate:
Unconjugated

Full name:
myosin, light chain 6, alkali, smooth muscle and non-muscle

Gene ID:
4637

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51759-1 targets MYL6 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MYL6 Monoclonal antibody, PBS Only (Capture) 68142-1-PBS (2G4C9). 100 μ g. Concentration 1 mg/mL.

Detection antibody: MYL6 Monoclonal antibody, PBS Only (Detector) 68142-2-PBS (2B3D6). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51759-1, MYL6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68142-1-PBS. Detection antibody: 68142-2-PBS.



Cytometric bead array standard curve of MP51759-1, MYL6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68142-1-PBS. Detection antibody: 68142-2-PBS. Standard: Ag13091. Range: 0.781-100 ng/mL.