

CEP55 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51760-1**

Capture Antibody Information

Catalog Number: 67460-2-PBS	Clone ID: 1A1B7	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: centrosomal protein 55kDa
Isotype: IgG1	GenBank: BC008947	Gene ID: 55165
Purification Method: Protein G purification	Immunogen Catalog Number: Ag29995	

Detection Antibody Information

Catalog Number: 67460-3-PBS	Clone ID: 2G9G4	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: centrosomal protein 55kDa
Isotype: IgG1	GenBank: BC008947	Gene ID: 55165
Purification Method: Protein G purification	Immunogen Catalog Number: Ag29995	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP51760-1 targets CEP55 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CEP55 Monoclonal antibody, PBS Only (Capture) 67460-2-PBS (1A1B7). 100 µg. Concentration 1 mg/mL.

Detection antibody: CEP55 Monoclonal antibody, PBS Only (Detector) 67460-3-PBS (2G9G4). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51760-1, CEP55 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67460-2-PBS. Detection antibody: 67460-3-PBS.



Cytometric bead array standard curve of MP51760-1, CEP55 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67460-2-PBS. Detection antibody: 67460-3-PBS. Standard: Ag29995. Range: 0.781-100 ng/mL.