

CHMP7 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51789-1**

Capture Antibody Information

Catalog Number: 68406-2-PBS	Clone ID: 2F1H6	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: CHMP family, member 7
Isotype: IgG1	GenBank: BC019110	Gene ID: 91782
Purification Method: Protein G purification	Immunogen Catalog Number: Ag9641	

Detection Antibody Information

Catalog Number: 68406-3-PBS	Clone ID: 1E8D5	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: CHMP family, member 7
Isotype: IgG2b	GenBank: BC019110	Gene ID: 91782
Purification Method: Protein G purification	Immunogen Catalog Number: Ag9641	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	--

Product Information

MP51789-1 targets CHMP7 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CHMP7 Monoclonal antibody, PBS Only (Capture) 68406-2-PBS (2F1H6). 100 μ g. Concentration 1 mg/mL.

Detection antibody: CHMP7 Monoclonal antibody, PBS Only (Detector) 68406-3-PBS (1E8D5). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51789-1, CHMP7 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68406-2-PBS. Detection antibody: 68406-3-PBS.



Cytometric bead array standard curve of MP51789-1, CHMP7 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68406-2-PBS. Detection antibody: 68406-3-PBS. Standard:Ag9641. Range: 0.781-100 ng/mL.