

CACYBP Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51795-2**

Capture Antibody Information

Catalog Number: 68161-1-PBS	Clone ID: 3C2G11	Conjugate: Unconjugated
Host: Mouse	Reactivity: human, mouse, rat	Full name: calcyclin binding protein
Isotype: IgG2a	GenBank: BC022352	Gene ID: 27101
Purification Method: Protein A purification	Immunogen Catalog Number: Ag2340	

Detection Antibody Information

Catalog Number: 68161-4-PBS	Clone ID: 3D5F2	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: calcyclin binding protein
Isotype: IgG2b	GenBank: BC022352	Gene ID: 27101
Purification Method: Protein A purification	Immunogen Catalog Number: Ag2340	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP51795-2 targets CACYBP in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CACYBP Monoclonal antibody, PBS Only (Capture) 68161-1-PBS (3C2G11). 100 µg. Concentration 1 mg/mL.

Detection antibody: CACYBP Monoclonal antibody, PBS Only (Detector) 68161-4-PBS (3D5F2). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data

