

EIF3H Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP51804-1

Capture Antibody Information

Catalog Number:
61110-1-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1B9G3

Reactivity:
human

GenBank:
BC000386

Immunogen Catalog Number:
Ag30831

Conjugate:
Unconjugated

Full name:
eukaryotic translation initiation factor 3, subunit H

Gene ID:
8667

Detection Antibody Information

Catalog Number:
61110-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
1G9D1

Reactivity:
human

GenBank:
BC000386

Immunogen Catalog Number:
Ag30831

Conjugate:
Unconjugated

Full name:
eukaryotic translation initiation factor 3, subunit H

Gene ID:
8667

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51804-1 targets EIF3H in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: EIF3H Monoclonal antibody, PBS Only (Capture) 61110-1-PBS (1B9G3). 100 μ g. Concentration 1 mg/mL.

Detection antibody: EIF3H Monoclonal antibody, PBS Only (Detector) 61110-2-PBS (1G9D1). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51804-1, EIF3H Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61110-1-PBS. Detection antibody: 61110-2-PBS.



Cytometric bead array standard curve of MP51804-1, EIF3H Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61110-1-PBS. Detection antibody: 61110-2-PBS. Standard:Ag30831. Range: 0.781-100 ng/mL