

EFTUD2 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51817-1**

Capture Antibody Information

Catalog Number:
67855-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2F4G2
Reactivity:
human, mouse, rat
GenBank:
BC002360
Immunogen Catalog Number:
Ag24345

Conjugate:
Unconjugated
Full name:
elongation factor Tu GTP binding domain containing 2
Gene ID:
9343

Detection Antibody Information

Catalog Number:
67855-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1H9G9
Reactivity:
human
GenBank:
BC002360
Immunogen Catalog Number:
Ag24345

Conjugate:
Unconjugated
Full name:
elongation factor Tu GTP binding domain containing 2
Gene ID:
9343

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51817-1 targets EFTUD2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: EFTUD2 Monoclonal antibody, PBS Only (Capture) 67855-1-PBS (2F4G2). 100 µg. Concentration 1 mg/mL.

Detection antibody: EFTUD2 Monoclonal antibody, PBS Only (Detector) 67855-2-PBS (1H9G9). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

