

FABP1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51819-1**

Capture Antibody Information

Catalog Number: 68227-2-PBS	Clone ID: 1A7D7	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: fatty acid binding protein 1, liver
Isotype: IgG1	GenBank: BC032801	Gene ID: 2168
Purification Method: Protein G purification	Immunogen Catalog Number: Ag32836	

Detection Antibody Information

Catalog Number: 68227-1-PBS	Clone ID: 1B6C8	Conjugate: Unconjugated
Host: Mouse	Reactivity: human, mouse, rat, pig	Full name: fatty acid binding protein 1, liver
Isotype: IgG1	GenBank: BC032801	Gene ID: 2168
Purification Method: Protein G purification	Immunogen Catalog Number: Ag32836	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP51819-1 targets FABP1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FABP1 Monoclonal antibody, PBS Only (Capture) 68227-2-PBS (1A7D7). 100 µg. Concentration 1 mg/mL.

Detection antibody: FABP1 Monoclonal antibody, PBS Only (Detector) 68227-1-PBS (1B6C8). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51819-1, FABP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68227-2-PBS. Detection antibody: 68227-1-PBS.



Cytometric bead array standard curve of MP51819-1, FABP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 6B227-2-PBS. Detection antibody: 6B227-1-PBS. Standard: Ag32836. Range: 0.781-100 ng/mL.