

SND1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP51825-1

Capture Antibody Information

Catalog Number:
60265-1-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
1A6A4

Reactivity:
human, mouse, rat

GenBank:
BC017180

Immunogen Catalog Number:
Ag1200

Conjugate:
Unconjugated

Full name:
staphylococcal nuclease and tudor domain containing 1

Gene ID:
27044

Detection Antibody Information

Catalog Number:
60265-2-PBS

Host:
Mouse

Isotype:
IgG2b

Purification Method:
Protein A purification

Clone ID:
2D5F6

Reactivity:
human

GenBank:
BC017180

Immunogen Catalog Number:
Ag1200

Conjugate:
Unconjugated

Full name:
staphylococcal nuclease and tudor domain containing 1

Gene ID:
27044

Applications

Tested Applications:
Cytometric bead array

Range:
0.156-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51825-1 targets SND1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SND1 Monoclonal antibody, PBS Only (Capture) 60265-1-PBS (1A6A4). 100 μ g. Concentration 1 mg/mL.

Detection antibody: SND1 Monoclonal antibody, PBS Only (Detector) 60265-2-PBS (2D5F6). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51825-1, SND1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60265-1-PBS. Detection antibody: 60265-2-PBS.



Cytometric bead array standard curve of MP51825-1, SND1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60265-1-PBS. Detection antibody: 60265-2-PBS. Standard:Ag1200. Range: 0.156-20 ng/mL