

Human IgG Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51853-1**

Capture Antibody Information

Catalog Number:
67761-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
4D11F7
Reactivity:
human
GenBank:

Conjugate:
Unconjugated
Full name:
Human IgG Kappa chain
Gene ID:

Detection Antibody Information

Catalog Number:
68275-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
1C1H2
Reactivity:
human
GenBank:
Immunogen Catalog Number:
Eg0212

Conjugate:
Unconjugated
Full name:
FC of Human IgG
Gene ID:

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51853-1 targets Human IgG in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Human IgG Kappa chain Monoclonal antibody, PBS Only (Capture) 67761-1-PBS (4D11F7). 100 μ g. Concentration 1 mg/mL.

Detection antibody: FC of Human IgG Monoclonal antibody, PBS Only (Detector) 68275-2-PBS (1C1H2). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51853-1, Human IgG Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67761-1-PBS. Detection antibody: 68275-2-PBS.



Cytometric bead array standard curve of MP51853-1, Human IgG Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67761-1-PBS. Detection antibody: 68275-2-PBS. Standard: Human IgG. Range: 0.098-100 ng/mL.