

BID Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51897-1**

Capture Antibody Information

Catalog Number:
60301-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein A purification

Clone ID:
2A2E5
Reactivity:
human
GenBank:
BC009197
Immunogen Catalog Number:
Ag21005

Conjugate:
Unconjugated
Full name:
BH3 interacting domain death agonist
Gene ID:
637

Detection Antibody Information

Catalog Number:
60301-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
4C7G2
Reactivity:
human
GenBank:
BC009197
Immunogen Catalog Number:
Ag21005

Conjugate:
Unconjugated
Full name:
BH3 interacting domain death agonist
Gene ID:
637

Applications

Tested Applications:
Cytometric bead array

Range:
0.391-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51897-1 targets BID in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: BID Monoclonal antibody, PBS Only (Capture) 60301-1-PBS (2A2E5). 100 µg. Concentration 1 mg/mL.

Detection antibody: BID Monoclonal antibody, PBS Only (Detector) 60301-2-PBS (4C7G2). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51897-1, BID Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60301-1-PBS. Detection antibody: 60301-2-PBS.



Cytometric bead array standard curve of MP51897-1, BID Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60301-1-PBS. Detection antibody: 60301-2-PBS. Standard:Ag21005. Range: 0.391-50 ng/mL.