

For Research Use Only

IL-1 beta Hybrid Matched Antibody Pair, PBS Only



Catalog Number:MP80001-4

Capture Antibody Information

Catalog Number:
82696-5-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
2E22
Reactivity:
Human
GenBank:
BC008678

Conjugate:
Unconjugated
Full name:
interleukin 1, beta
Gene ID:
3553

Detection Antibody Information

Catalog Number:
68767-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
3H8A3
Reactivity:
Human
GenBank:
BC008678

Conjugate:
Unconjugated
Full name:
interleukin 1, beta
Gene ID:
3553

Applications

Tested Applications:
Cytometric bead array

Range:
3.125-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP80001-4 targets IL-1 beta in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: IL1B Recombinant antibody, PBS Only (Capture/Detector) 82696-5-PBS (2E22). 100 μ g. Concentration 1 mg/mL.

Detection antibody: IL-1 beta Monoclonal antibody, PBS Only (Detector) 68767-2-PBS (3H8A3). 100 μ g. Concentration 1 mg/mL.

Alternative IL-1 beta matched antibody pairs: MP00331-1, MP50134-1, MP80001-1, MP80001-2, MP80001-3, MP80001-5, MP80001-6

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP80001-4, IL-1 beta Hybrid Matched Antibody Pair, PBS Only. Capture antibody: 82696-5-RR. Detection antibody: 68767-2-PBS. Standard:HZ-1164. Range:3.125-100 ng/mL