

For Research Use Only

FcZero-rAb® PE-CoraLite® Plus 594 Anti-Human CD2 (TS1/8) Rabbit Recombinant Antibody

Catalog Number: PCL594-FcA65565



Basic Information

Catalog Number:
PCL594-FcA65565

Concentration:
100tests, 5 ul/test

Source:
Rabbit

Isotype:
IgG

GenBank Accession Number:
BC033583

GeneID (NCBI):
914

ENSEMBL Gene ID:
ENSG00000116824

UNIPROT ID:
P06729

Full Name:
CD2 molecule

Calculated MW:
351 aa, 39 kDa

Purification Method:
Protein A purification

CloneNo.:
TS1/8

Recommended Dilutions:
FC: 5 ul per 10⁶ cells in a 100 µl suspension

Excitation/Emission maxima wavelengths:
496 nm, 565 nm / 604 nm

Applications

Tested Applications:
FC

Species Specificity:
human

Positive Controls:

FC : human PBMCs,

Background Information

CD2 is a cell surface glycoprotein present on a majority of thymocytes, all mature T cells and subset of NK cells but not on B lymphocytes. It is a pan T-cell marker. CD2 interacts with lymphocyte function-associated antigen (LFA-3/CD58) and CD48/BCM1 to mediate adhesion between T-cells and other cell types. CD2 is implicated in the triggering of T-cells, the cytoplasmic domain is implicated in the signaling function.

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ human PBMCs were surface stained with FcZero-rAb® CoraLite® Plus 488 Anti-Human CD19 (4G7), and 5 ul PE-CoraLite® Plus 594 Anti-Human CD2 (TS1/8) Rabbit RecAb (PCL594-FcA65565, Clone: TS1/8) or FcZero-rAb® PE-CoraLite® Plus 594 Rabbit IgG Isotype Control Recombinant Antibody (PCL594-FcA98136, Clone: 240953C9). Cells were not fixed. Lymphocytes were gated.



1x10⁶ human PBMCs were surface stained with 5 ul PE-CoraLite® Plus 594 Anti-Human CD2 (TS1/8) Rabbit RecAb (PCL594-FcA65565, Clone: TS1/8) (red) or FcZero-rAb® PE-CoraLite® Plus 594 Rabbit IgG Isotype Control Recombinant Antibody (PCL594-FcA98136, Clone: 240953C9) (blue). Cells were not fixed. Lymphocytes were gated.