

For Research Use Only

FcZero-rAb™ PE-CoraLite® Plus 594 Anti-Human CD226 (11A8) Rabbit Recombinant Antibody

Catalog Number: PCL594-FcA65573



Basic Information

Catalog Number:
PCL594-FcA65573

Concentration:
100 tests, 5 µl/test

Source:
Rabbit

Isotype:
IgG

GenBank Accession Number:
BC074787

GeneID (NCBI):
10666

ENSEMBL Gene ID:
ENSG00000150637

UNIPROT ID:
Q15762

Full Name:
CD226 molecule

Calculated MW:
336 aa, 39 kDa

Purification Method:
Protein A purification

CloneNo.:
11A8

Recommended Dilutions:
FC: 5 µl per 10⁶ cells in a 100 µl suspension

Excitation/Emission maxima wavelengths:
496 nm, 565 nm / 604 nm

Applications

Tested Applications:
FC

Species Specificity:
human

Positive Controls:

FC : human PBMCs,

Background Information

CD226 (DNAM-1) is a ~65 kDa glycoprotein expressed on the surface of NK cells, platelets, monocytes and a subset of T cells. It is a member of the Ig-superfamily containing 2 Ig-like domains of the V-set. CD226 mediates cellular adhesion of platelets and megakaryocytic cells to vascular endothelial cells. The protein also plays a role in megakaryocytic cell maturation. Interactions of CD226 and its ligands, CD155 and CD112, induce NK and T cell-mediated cytotoxicity and cytokine secretion (PMID: 15039383).

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ human PBMCs were surface stained with CoraLite® Plus 488 Anti-Human CD3 (UCHT1), and 5 ul PE-CoraLite® Plus 594 Anti-Human CD226 (11A8) Rabbit IgG RecAb (PCL594-FcA65573, Clone: 11A8) or Rabbit IgG Isotype Control RecAb (PCL594-FcA98136, Clone: 240953C9). Cells were not fixed.



1x10⁶ human PBMCs were surface stained with 5 ul PE-CoraLite® Plus 594 Anti-Human CD226 (11A8) Rabbit IgG RecAb (PCL594-FcA65573, Clone: 11A8) (red) or Rabbit IgG Isotype Control RecAb (PCL594-FcA98136, Clone: 240953C9) (blue). Cells were not fixed.