

For Research Use Only

PE-conjugated TXNDC5 Polyclonal antibody

Catalog Number: PE-19834

Featured Product



Basic Information

Catalog Number:

PE-19834

Concentration:

200 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13808

GenBank Accession Number:

BC001199

GeneID (NCBI):

81567

UNIPROT ID:

Q8NBS9

Full Name:

thioredoxin domain containing 5
(endoplasmic reticulum)

Calculated MW:

432 aa, 48 kDa

Observed MW:

48 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IF/ICC: 1:50-1:500

Excitation/Emission maxima
wavelengths:

496 nm, 565 nm / 578 nm

Applications

Tested Applications:

IF/ICC

Species Specificity:

human, mouse, rat

Positive Controls:

IF/ICC : U2OS cells,

Background Information

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed U2OS cells using PE TXNDC5 antibody (PE-19834) at dilution of 1:200.