

PE Anti-Human PD-1/CD279 (J110)

Catalog Number: PE-65119

Basic Information

Catalog Number:

PE-65119

Size:

100 tests, 5 µl/test

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC074740

GeneID (NCBI):

5133

UNIPROT ID:

Q15116

Full Name:

programmed cell death 1

Calculated MW:

288 aa, 32 kDa

Purification Method:

Affinity purification

CloneNo.:

J110

Excitation/Emission maxima wavelengths:

496 nm, 565 nm / 578 nm

Applications

Tested Applications:

FC

Species Specificity:

human

Background Information

Programmed cell death 1 (PD-1, also known as CD279) is an immunoinhibitory receptor that belongs to the CD28/CTLA-4 subfamily of the Ig superfamily. It is a 288 amino acid (aa) type I transmembrane protein composed of one Ig superfamily domain, a stalk, a transmembrane domain, and an intracellular domain containing an immunoreceptor tyrosine-based inhibitory motif (ITIM) as well as an immunoreceptor tyrosine-based switch motif (ITSM) (PMID: 18173375). PD-1 is expressed during thymic development and is induced in a variety of hematopoietic cells in the periphery by antigen receptor signaling and cytokines (PMID: 20636820). Engagement of PD-1 by its ligands PD-L1 or PD-L2 transduces a signal that inhibits T-cell proliferation, cytokine production, and cytolytic function (PMID: 19426218). It is critical for the regulation of T cell function during immunity and tolerance. Blockade of PD-1 can overcome immune resistance and also has been shown to have antitumor activity (PMID: 22658127; 23169436).

Storage


Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.


Storage Buffer:

Phosphate based buffer with 0.09% sodium azide and 0.1% gelatin, pH 7.2.

Selected Validation Data



1x10⁶ untreated or PHA-treated human PBMCs were surface stained with APC Anti-Human CD3 and 5 ul PE Anti-Human PD-1/CD279 (J110) (PE-65119, Clone: J110). Cells were not fixed.



1x10⁶ PHA treated human PBMCs were surface stained with 5 ul PE Anti-Human PD-1/CD279 (J110) (PE-65119, Clone: J110) (red) or PE Mouse IgG1 Isotype Control (MOPC-21) (PE-65124, Clone: MOPC-21) (blue). 1x10⁶ untreated human PBMCs were surface stained with 5 ul PE Anti-Human PD-1/CD279 (J110) (PE-65119, Clone: J110) (black, dashed). Cells were not fixed.