

For Research Use Only

PE Anti-Human CD120a (H398)

Catalog Number: **PE-65170**



Basic Information

Catalog Number:

PE-65170

Size:

100 tests, 10 µl/test

Source:

Mouse

Isotype:

IgG2a

GenBank Accession Number:

BC010140

GeneID (NCBI):

7132

UNIPROT ID:

P19438

Full Name:

tumor necrosis factor receptor
superfamily, member 1A

Calculated MW:

455 aa, 50 kDa

Purification Method:

The purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. The conjugate is purified by size-exclusion chromatography.

CloneNo.:

H398

Excitation/Emission maxima
wavelengths:

496 nm, 565 nm / 578 nm

Applications

Tested Applications:

FC

Species Specificity:

Human

Background Information

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide.

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



100 ul human peripheral blood were surface stained with 10.00 ul PE Anti-Human CD120a (PE-65170, Clone: H398). Cells were not fixed.



100 ul human peripheral blood were surface stained with PE-mouse IgG2a isotype control antibody. Cells were not fixed.