

For Research Use Only

PE Anti-Human CD1c (L161)

Catalog Number: **PE-65176**



Basic Information

Catalog Number:

PE-65176

Size:

100 tests, 10 µl/test

Source:

Mouse

Isotype:

Mouse IgG1 kappa

GenBank Accession Number:

BC126465

GeneID (NCBI):

911

UNIPROT ID:

P29017

Full Name:

CD1c molecule

Purification Method:

The purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. The conjugate is purified by size-exclusion chromatography.

CloneNo.:

L161

Excitation/Emission maxima wavelengths:

496 nm, 565 nm / 578 nm

Applications

Tested Applications:

FC

Species Specificity:

Human

Background Information

Storage

Storage:

Store at 2-8°C. Avoid exposure to light.

Storage Buffer:

PBS with 0.1% sodium azide.

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ human peripheral blood lymphocytes were surface stained with APC Anti-Human CD19 (APC-65110, Clone:HI B19) and 10.00 ul/test PE-Anti-Human CD1c (PE-65176, Clone:L161). Cells were not fixed.



1X10⁶ human peripheral blood lymphocytes were surface stained with APC Anti-Human CD19 (APC-65110, Clone:HI B19) and PE-mouse IgG1 isotype control antibody. Cells were not fixed.