

For Research Use Only

# PE Anti-Human CD81 (5A6)

Catalog Number: PE-65195



## Basic Information

Catalog Number:

PE-65195

Concentration:

100 tests, 5 µl/test

Source:

Mouse

Isotype:

IgG1, kappa

GenBank Accession Number:

BC002978

GeneID (NCBI):

975

ENSEMBL Gene ID:

ENSG00000110651

UNIPROT ID:

P60033

Full Name:

CD81 molecule

Calculated MW:

26 kDa

Purification Method:

Affinity purification

CloneNo.:

5A6

Excitation/Emission maxima wavelengths:

496 nm, 565 nm / 578 nm

## Applications

Tested Applications:

FC

Species Specificity:

human

## Background Information

CD81 (also known as TAPA1 or TSPAN28) is a membrane protein of the tetraspanin superfamily, which are characterized by the presence of four conserved transmembrane regions. Many of these members are expressed on leukocytes and have been implicated in signal transduction, cell-cell interactions, and cellular activation and development. CD81 is involved in signal transduction and cell adhesion in the immune system (PMID: 9597125). CD81 has also been identified as an essential receptor for HCV (hepatitis C virus) (PMID: 21428934).

## Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1x10<sup>6</sup> human PBMCs were surface stained with 5 ul PE Anti-Human CD81 (PE-65195, Clone:5A6) or PE Mouse IgG1 Isotype Control (PE-65124). Cells were not fixed. Lymphocytes were gated.



1x10<sup>6</sup> human PBMCs were surface stained with 5 ul PE Anti-Human CD81 (PE-65195, Clone:5A6) or PE Mouse IgG1 Isotype Control (PE-65124). Cells were not fixed. Lymphocytes were gated.