

For Research Use Only

# PE Anti-Human CD11c (BU15)

Catalog Number: PE-65196



## Basic Information

**Catalog Number:**

PE-65196

**Size:**

100 tests, 5  $\mu$  l / test

**Source:**

Mouse

**Isotype:**

IgG1, kappa

**GenBank Accession Number:**

BC038237

**GeneID (NCBI):**

3687

**ENSEMBL Gene ID:**

ENSG00000140678

**UNIPROT ID:**

P20702

**Full Name:**

integrin, alpha X (complement component 3 receptor 4 subunit)

**Calculated MW:**

1169 aa, 129 kDa

**Purification Method:**

Affinity purification

**CloneNo.:**

BU15

**Excitation/Emission maxima wavelengths:**

496 nm, 565 nm / 578 nm

## Applications

**Tested Applications:**

FC

**Species Specificity:**

Human

## Background Information

Integrins are cell adhesion receptors that are heterodimers composed of non-covalently associated  $\alpha$  and  $\beta$  subunits (PMID: 9779984). CD11c, also known as integrin  $\alpha$  X, is a type I transmembrane glycoprotein present on a variety of cells, including monocytes/macrophages, granulocytes, a subset of B cells, NK cells and dendritic cells (PMID: 2897326; 1680915; 1694698; 17389580). As a result of its high level of expression on most dendritic cells, CD11c is typically considered to be a marker of conventional dendritic cells (PMID: 27119555). CD11c forms an  $\alpha$  /  $\beta$  heterodimer with CD18 (integrin  $\beta$  2). CD11c/CD18 acts as a receptor for fibrinogen and is important in monocyte adhesion and chemotaxis (PMID: 1671533).

## Storage

**Storage:**

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

**Storage Buffer:**

PBS with 0.09% sodium azide and 0.5% BSA.

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> human PBMCs were incubated with FC receptor block then surface stained with 5 ul PE Anti-Human CD11c (PE-65196, Clone:BU15), or PE Mouse IgG1 Isotype Control. Cells were not fixed. Monocytes were gated.



1X10<sup>6</sup> human PBMCs were incubated with FC receptor block then surface stained with 5 ul PE Anti-Human CD11c (PE-65196, Clone:BU15), or PE Mouse IgG1 Isotype Control. Cells were not fixed. Monocytes were gated.