

For Research Use Only

# FcZero-rAb® PE Anti-Mouse Ly-6C (HK1.4) Rabbit Recombinant Antibody

Catalog Number: PE-FcA65669



## Basic Information

Catalog Number:

PE-FcA65669

Concentration:

100ug, 200 ug/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC010764

GeneID (NCBI):

17067

UNIPROT ID:

POCW02

Full Name:

lymphocyte antigen 6 complex, locus C1

Purification Method:

Protein A purification

CloneNo.:

HK1.4

Recommended Dilutions:

FC: 0.1 ug per 10<sup>6</sup> cells in a 100 µl suspension

Excitation/Emission maxima wavelengths:

496 nm, 565 nm / 578 nm

## Applications

Tested Applications:

FC

Species Specificity:

mouse

Positive Controls:

FC : mouse bone marrow cells,

## Background Information

Ly-6C (lymphocyte antigen 6 complex, locus C) is a 14-17 kDa GPI-linked cell surface antigen that belongs to the Ly-6 family of glycosyl phosphatidylinositol-anchored surface glycoproteins (PMID: 9192663). Ly-6C is expressed on mouse monocytes/macrophages, endothelial cells, granulocytes, thymocytes, and some T cell subsets (PMID: 2849552). The two alleles of Ly-6C, Ly-6C.1 and Ly-6C.2, result in differential expression of this protein on subpopulations of CD4<sup>+</sup> T cells (PMID: 15057973; 9126982). Ly-6C has been reported as a marker of memory CD8<sup>+</sup> T cells (PMID: 7543536).

## Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1x10<sup>6</sup> mouse bone marrow cells were surface stained with FcZero-rAb® Coralite® Plus 647 Anti-Mouse CD11b, and 0.1 ug PE Anti-Mouse Ly-6C (HK1.4) Rabbit RecAb (PE-FcA65669, Clone: HK1.4) or FcZero-rAb® PE Rabbit IgG Isotype Control Recombinant Antibody (PE-FcA98136, Clone: 240953C9). Cells were not fixed.



1x10<sup>6</sup> mouse bone marrow cells were surface stained with 0.1 ug PE Anti-Mouse Ly-6C (HK1.4) Rabbit RecAb (PE-FcA65669, Clone: HK1.4) (red) or FcZero-rAb® PE Rabbit IgG Isotype Control Recombinant Antibody (PE-FcA98136, Clone: 240953C9) (blue). Cells were not fixed.