

For Research Use Only

# FcZero-rAb® PE Anti-Human CD34 Rabbit Recombinant Antibody

Catalog Number: PE-FcA98145-4



## Basic Information

### Catalog Number:

PE-FcA98145-4

### Concentration:

100 tests, 5 µl/test

### Source:

Rabbit

### Isotype:

IgG

### GenBank Accession Number:

BC039146

### GeneID (NCBI):

947

### UNIPROT ID:

P28906

### Full Name:

CD34 molecule

### Calculated MW:

41 kDa

### Purification Method:

Protein A purification

### CloneNo.:

242069A6

### Recommended Dilutions:

FC: 5 µl per 10<sup>6</sup> cells in a 100 µl suspension

### Excitation/Emission maximum wavelengths:

496 nm, 565 nm / 578 nm

## Applications

### Tested Applications:

FC

### Species Specificity:

human

### Positive Controls:

FC : Transfected HEK-293T cells,

## Background Information

CD34 is a 105- to 120-kDa glycoprophosphoprotein expressed on the majority of hematopoietic stem/progenitor cells, bone marrow stromal cells, capillary endothelial cells, embryonic fibroblasts, and some nerve tissue. CD34 is a commonly used marker for identifying human hematopoietic stem/progenitor cells and mediates cell adhesion and lymphocyte homing by binding L-selectin and E-selectin ligands. CD34 is also one of the best negative selection markers for characterizing and/or isolating human MSCs from bone marrow and other sources. Along with other positive selection markers (such as CD29, CD44, CD90, CD105 and CD166), negative selection markers (such as CD34 and CD45) are used for MSC identification.

## Storage

### Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

### Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1x10<sup>6</sup> HEK-293T cells (left) or CD34-transfected HEK-293T cells (right) were surface stained with 5 ul PE Anti-Human CD34 Rabbit RecAb (PE-FcA98145-4, Clone: 242069A6). Cells were not fixed.



1x10<sup>6</sup> HEK-293T cells (blue) or CD34-transfected HEK-293T cells (red) were surface stained with 5 ul PE Anti-Human CD34 Rabbit RecAb (PE-FcA98145-4, Clone: 242069A6). Cells were not fixed.