

For Research Use Only

FcZero-rAb™ PE Anti-Mouse LAG-3/CD223 Rabbit Recombinant Antibody

Catalog Number: PE-FcA98167



Basic Information

Catalog Number:

PE-FcA98167

Concentration:

100ug, 200 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG1220

GenBank Accession Number:

NM_008479.2

GeneID (NCBI):

16768

Full Name:

Lymphocyte-activation gene 3

Calculated MW:

57 kDa

Purification Method:

Protein A purification

CloneNo.:

241294E11

Recommended Dilutions:

FC: 0.1 ug per 10⁶ cells in a 100 µl suspension

Excitation/Emission maximum wavelengths:

496 nm, 565 nm / 578 nm

Applications

Tested Applications:

FC

Species Specificity:

mouse

Positive Controls:

FC : Anti-CD3/CD28 treated mouse splenocytes,

Background Information

Lymphocyte activation gene 3 protein (LAG-3), also known as CD223, is a type I transmembrane protein belonging to the immunoglobulin superfamily (PMID: 1692078). It is a CD4-related molecule that binds MHC class II with high affinity (PMID: 1692078; 9159144). LAG-3 is expressed on activated T and NK cells. Besides, natural CD4(+)CD25(+) Tregs express LAG-3 upon activation (PMID: 15485628). LAG-3 negatively regulates T-cell function and homeostasis, and contributes to the suppressor activity of Tregs (PMID: 14644131; 15485628).

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ untreated or anti-CD3/CD28 treated mouse splenocytes were surface stained with FcZero-rAb™ APC Anti-Mouse CD3e and 0.1 ug PE Anti-Mouse LAG-3/CD223 Rabbit RecAb (PE-FcA98167, Clone: 241294E11). Cells were not fixed.