

For Research Use Only



FcZero-rAb® PE Anti-Human DC-SIGNR/CD299 Rabbit Recombinant Antibody

Catalog Number: PE-FcA98205

Basic Information

Catalog Number: PE-FcA98205	GenBank Accession Number: BC038851	Purification Method: Protein A purification
Concentration: 100 tests, 5 µl/test	GeneID (NCBI): 10332	CloneNo.: 241683G8
Source: Rabbit	UNIPROT ID: Q9H2X3	Recommended Dilutions: FC: 5 µl per 10 ⁶ cells in a 100 µl suspension
Isotype: IgG	Full Name: C-type lectin domain family 4, member M	Excitation/Emission maxima wavelengths: 496 nm, 565 nm / 578 nm
	Calculated MW: 45 kDa	

Applications

Tested Applications: FC	Positive Controls: FC : Transfected HEK-293T cells,
Species Specificity: human	

Background Information

Dendritic cell-specific intracellular adhesion molecule-3-grabbing non-integrin (DC-SIGN/CD209) and DC-SIGN related (DC-SIGNR/CD299/L-SIGN/CD209L/CLEC4M) are trans-membrane C-type lectins that are calcium-dependent carbohydrate-binding proteins, which can act as cell-adhesion and pathogen-recognition receptors (PMID: 18458800). DC-SIGN and DC-SIGNR share 77% amino acid identity (PMID: 10975799; 11226297). DC-SIGN is expressed on subsets of dendritic cells and macrophages while DC-SIGNR is restricted to endothelial cells of the liver, lymph nodes and placenta (PMID: 17509452). DC-SIGNR serves as a receptor for ICAM-3. It recognizes a wide range of evolutionarily divergent pathogens.

Storage

Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 0.09% sodium azide and 0.5% BSA, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ HEK-293T cells or DC-SIGNR/CD299 transfected HEK-293T cells were surface stained with 5 μ l PE Anti-Human DC-SIGNR/CD299 Rabbit RecAb (PE-FcA98205, Clone: 241683G8). Cells were not fixed.