

For Research Use Only

# FcZero-rAb® PE-Cyanine5 Anti-Human CD3 (OKT3) Rabbit Recombinant Antibody

Catalog Number: PY5-FcA65569



## Basic Information

Catalog Number:

PY5-FcA65569

Concentration:

100tests, 5 ul/test

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC049847

GeneID (NCBI):

916

ENSEMBL Gene ID:

ENSG00000198851

UNIPROT ID:

P07766

Full Name:

CD3e molecule, epsilon (CD3-TCR complex)

Calculated MW:

207 aa, 23 kDa

Purification Method:

Protein A purification

CloneNo.:

OKT3

Recommended Dilutions:

FC: 5 ul per 10<sup>6</sup> cells in a 100 µl suspension

Excitation/Emission maximum wavelengths:

496 nm, 565 nm / 666 nm

## Applications

Tested Applications:

FC

Species Specificity:

human

Positive Controls:

FC : human PBMCs,

## Background Information

CD3 is a multimeric protein associated with the T-cell receptor (TCR) to form a complex involved in antigen recognition and signal transduction (PMID: 15885124). CD3 is composed of CD3  $\gamma$ ,  $\delta$ ,  $\epsilon$ , and  $\zeta$  chains (PMID: 1826255). It is expressed by thymocytes in a developmentally regulated manner, T cells, and some NK cells (PMID: 3289580). The TCR recognizes antigens bound to major histocompatibility complex (MHC) molecules. TCR-mediated peptide-MHC recognition is transmitted to the CD3 complex, leading to the intracellular signal transduction (PMID: 11985657).

## Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1x10<sup>6</sup> human PBMCs were surface stained with FcZero-rAb™ CoraLite® Plus 488 Anti-Human CD19 (4G7), and 5 ul PE-Cyanine5 Anti-Human CD3 (OKT3) Rabbit RecAb (PY5-FcA65569, Clone: OKT3) or FcZero-rAb™ PE-Cyanine5 Rabbit IgG Isotype Control Recombinant Antibody (PY5-FcA98136, Clone: 240953C9). Cells were not fixed. Lymphocytes were gated.



1x10<sup>6</sup> human PBMCs were surface stained with 5 ul PE-Cyanine5 Anti-Human CD3 (OKT3) Rabbit RecAb (PY5-FcA65569, Clone: OKT3) (red) or FcZero-rAb™ PE-Cyanine5 Rabbit IgG Isotype Control Recombinant Antibody (PY5-FcA98136, Clone: 240953C9) (blue). Cells were not fixed. Lymphocytes were gated.