

For Research Use Only

# FcZero-rAb® PE-Cyanine5 Anti-Rat CD4 Rabbit Recombinant Antibody

Catalog Number: PY5-FcA98252



## Basic Information

Catalog Number:

PY5-FcA98252

Concentration:

100ug, 200 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2025

GenBank Accession Number:

NM\_012705.1

GeneID (NCBI):

24932

UNIPROT ID:

P05540

Full Name:

Cd4 molecule

Calculated MW:

51 kDa

Purification Method:

Protein A purification

CloneNo.:

241955C11

Recommended Dilutions:

FC: 0.1 ug per 10<sup>6</sup> cells in a 100 µl suspension

Excitation/Emission maximum wavelengths:

496 nm, 565 nm / 666 nm

## Applications

Tested Applications:

FC

Species Specificity:

rat

Positive Controls:

FC : rat splenocytes,

## Background Information

CD4 is a 55-kDa transmembrane glycoprotein expressed on T helper cells, majority of thymocytes, monocytes, macrophages, and dendritic cells (PMID: 9304802; 12213222). CD4 is an accessory protein for MHC class-II antigen/T-cell receptor interaction. It plays an important role in T helper cell development and activation (PMID: 9539765; 3112582). CD4 serves as a receptor for the human immunodeficiency virus (HIV) (PMID: 9304802).

## Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1x10<sup>6</sup> rat splenocytes were surface stained with 0.1 ug PE-Cyanine5 Anti-Rat CD4 Rabbit RecAb (PY5-FcA98252, Clone: 241955C11) or FcZero-rAb™ PE-Cyanine5 Rabbit IgG Isotype Control Recombinant Antibody (PY5-FcA98136, Clone: 240953C9). Cells were not fixed.



1x10<sup>6</sup> rat splenocytes were surface stained with 0.1 ug PE-Cyanine5 Anti-Rat CD4 Rabbit RecAb (PY5-FcA98252, Clone: 241955C11) (red) or FcZero-rAb™ PE-Cyanine5 Rabbit IgG Isotype Control Recombinant Antibody (PY5-FcA98136, Clone: 240953C9) (blue). Cells were not fixed.