

Multi-rAb™ Cleaved Caspase 3/P17/P19 Multi-Recombinant antibody

Catalog Number: **RMX00003**

Basic Information

Catalog Number: RMX00003	GenBank Accession Number: NM_032991	Purification Method: N/A
Concentration: 500 ug/ml	GeneID (NCBI): 836	Recommended Dilutions: WB 1:1000-1:6000 IHC 1:250-1:1000 IF/ICC 1:200-1:800
Source: Rabbit	UNIPROT ID: P42574	
Isotype: IgG	Full Name: caspase 3, apoptosis-related cysteine peptidase	
	Calculated MW: 32 kDa	
	Observed MW: 17 kDa, 19 kDa	

Applications

Tested Applications:
WB, IHC, IF/ICC, ELISASpecies Specificity:
human, mouse

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB : Jurkat cells,

IHC : mouse brain tissue,

IF/ICC : HeLa cells,

Background Information

Caspases, a family of endoproteases, are critical players in cell regulatory networks controlling inflammation and cell death. Initiator caspases (caspase-2, -8, -9, -10, -11, and -12) cleave and activate downstream effector caspases (caspase-3, -6, and -7), which in turn execute apoptosis by cleaving targeted cellular proteins. Caspase 3 plays a key role in the activation of sterol regulatory element binding proteins (SREBPs) between the basic helix-loop-helix leucine zipper domain and the membrane attachment domain. Caspase 3 can form heterocomplex with other proteins. This antibody can recognize p17 of Caspase 3. This antibody is specific for cleaved caspase 3, and does not recognize full length caspase-3, it can recognize the 17 kDa cleaved-caspase 3. The cleaved P17 and P19 fragment might form complex and shows at around 30-35 kDa by western blot (PMID: 25501826).

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Staurosporine treated and untreated Jurkat cells were subjected to SDS PAGE followed by western blot with RMX00003 (Cleaved Caspase 3/P17/P19) at dilution of 1:3000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using RMX00003 (Cleaved Caspase 3/P17/P19 antibody) at dilution of 1:500 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using RMX00003 (Cleaved Caspase 3/P17/P19 antibody) at dilution of 1:500 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using Cleaved Caspase 3/P17/P19 antibody (RMX00003) at dilution of 1:400 and Multi-rAb CoraLite ® Plus 488-Goat Anti-Rabbit Recombinant Secondary Antibody (H+L) (RGAR002).