

For Research Use Only.
Not For Use In Diagnostics.

HRP-conjugated Donkey Anti-Rat IgG(H+L)

Catalog Number: SA00001-10



Information

Catalog Number:

SA00001-10

Host:

Donkey

Applications:

ELISA, WB, Dot blot

Reactivity:

Rat (Minimal Cross Reactivity: Bovine, Chicken, Goat, Guinea Pig, Syrian Hamster, Horse, Human, Mouse, Rabbit, Sheep Serum Proteins)

Physical State:

Liquid

Conjugation:

HRP

Recommended Dilutions

1:2000-1:10,000 for western blotting with ECL substrates

1:1000-1:20,000 for ELISA and Western blotting with chromogenic substrates

Stabilizer

3 mg/ml BSA

Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 3 mg/ml BSA, pH 7.6.

Aliquoting is unnecessary for -20°C storage

Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.


For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Rat plasma and rat serum were subjected to SDS PAGE followed by western blot with SA00001-10 (HRP-conjugated Donkey Anti-Rat IgG(H+L)) at dilution of 1:10000 incubated at room temperature for 1.5 hours.