

For Research Use Only.  
Not For Use In Diagnostics.

# HRP -Protein A

Catalog Number: SA00001-18



[www.ptglab.com](http://www.ptglab.com)

## Information

Catalog Number:

SA00001-18

Host:

Applications:

WB, Dot blot

Reactivity:

FC Fragment

Physical State:

Liquid

Conjugation:

HRP

## Recommended Dilutions

1:4000-1:10000 for Western Blotting with chromogenic substrates after immunoprecipitation

## Stabilizer

3 mg/ml BSA

## Safety Notes

1.This product is for research use only,not for diagnostic or therapeutic use;<br>2. Do not use Sodium Azide as preservative.

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 3 mg/ml BSA, pH 7.6.

Aliquoting is unnecessary for -20°C storage

## Purity

Protein A was purified by chromatography from Staphylococcus aureus

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



IP Result of anti-EBAG9 (IP:12255-1-AP, 3ug;  
Detection:12255-1-AP 1:800) with mouse heart  
tissue lysate. SA00001-18 (HRP-conjugated Protein  
A from *S.aureus*) as secondary antibody.