

For Research Use Only.  
Not For Use In Diagnostics.

# HRP-conjugated Rabbit Anti-Goat IgG(H+L)

Catalog Number:SA00001-4



## Information

Catalog Number:  
SA00001-4  
Host:  
Rabbit  
Applications:  
ELISA, WB, Dot blot

Reactivity:  
Goat  
Physical State:  
Liquid  
Conjugation:  
HRP

## Recommended Dilutions

1:2000-1:10,000 for western blotting with ECL substrates  
1:1000-1:20,000 for ELISA and Western blotting with chromogenic substrates

## Stabilizer

3 mg/ml BSA

## Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 3 mg/ml BSA, pH 7.6.  
Aliquoting is unnecessary for -20°C storage


## Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please contact:  
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



IgG proteins of various species were subjected to SDS PAGE followed by western blot with SA00001-4 (HRP-conjugated Rabbit Anti-Goat IgG(H+L)) at dilution of 1:1000 incubated at room temperature for 1.5 hours.