

For Research Use Only.
Not For Use In Diagnostics.

HRP-conjugated Rabbit Anti-Chicken IgY (IgG)(H+L)

Catalog Number:SA00001-6



Information

Catalog Number:
SA00001-6
Host:
Rabbit
Applications:
ELISA, WB, Dot blot

Reactivity:
Chicken
Physical State:
Liquid
Conjugation:
HRP

Recommended Dilutions

1:2000-1:10,000 for western blotting with ECL substrates
1:1000-1:20,000 for ELISA and Western blotting with chromogenic substrates

Stabilizer

3 mg/ml BSA

Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 3 mg/ml BSA, pH 7.6.
Aliquoting is unnecessary

Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please contact:

T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

