

For Research Use Only.
Not For Use In Diagnostics.

Fluorescein (FITC)–conjugated Goat Anti-Human IgG(H+L)

Catalog Number:SA00003-12



Information

Catalog Number:

SA00003-12

Host:

Goat

Applications:

IF, FC, WB

Reactivity:

Human

Physical State:

Liquid

Conjugation:

FITC

Recommended Dilutions

1:20-1:100 for IF or Flow cyt; 1:1000-1:5000 for WB

Preservative

0.02% Sodium Azide

Fluorophore

Fluorescein-5-isothiocyanate (FITC) Amax=492nm; Emax=520nm

Fluorophore Protein

8.1 ug/mg; 3.1 moles FITC per mole IgG

Stabilizer

5 mg/ml BSA

Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 5 mg/ml BSA, 0.02% Sodium Azide, pH 7.6.

Aliquoting is unnecessary for -20°C storage

Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Human IgG (0.2 ug) were subjected to SDS PAGE followed by western blot with SA00003-12 at dilution of 1:1000 incubated at room temperature for 1.5 hours.