

For Research Use Only.  
Not For Use In Diagnostics.

# Fluorescein (FITC)-conjugated Goat Anti-Rabbit IgG(H+L)

Catalog Number:SA00003-2



[www.ptglab.com](http://www.ptglab.com)

## Information

Catalog Number:

SA00003-2

Host:

Goat

Applications:

IF, FC

Reactivity:

Rabbit

Physical State:

Liquid

Conjugation:

FITC

## Recommended Dilutions

1:100-1:1000 for most applications

## Preservative

0.02% Sodium Azide

## Fluorophore

Fluorescein-5-isothiocyanate (FITC)  $\lambda_{max}=492nm$ ;  $\lambda_{em}=520nm$

## Fluorophore Protein

9.0 ug/mg; 3.4 moles FITC per mole IgG

## Stabilizer

5 mg/ml BSA

## Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

## Storage

**Storage:**

Store at  $-20^{\circ}C$ . Stable for one year after shipment.

**Storage Buffer:**

0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 5 mg/ml BSA, 0.02% Sodium Azide, pH 7.6.

Aliquoting is unnecessary for  $-20^{\circ}C$  storage

## Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.


For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com) W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



1X10<sup>6</sup> HepG2 cells were intracellularly stained with 0.2 ug Anti-Human MYH9 (11128-1-AP) and Fluorescein (FITC)-conjugated Goat Anti-Rabbit IgG(H+L) (SA00003-2) at dilution 1:200 (red), and 0.2 ug Control Antibody. Cells were fixed with 90% MeOH.