

For Research Use Only.  
Not For Use In Diagnostics.

# Fluorescein (FITC)–conjugated Donkey Anti-Rabbit IgG(H+L)

Catalog Number:SA00003-8



## Information

### Catalog Number:

SA00003-8

### Host:

Donkey

### Applications:

IF, FC

### Reactivity:

Rabbit

### Physical State:

Liquid

### Conjugation:

FITC

## Recommended Dilutions

1:20-1:100 for most applications

## Preservative

0.02% Sodium Azide

## Fluorophore

Fluorescein-5-isothiocyanate (FITC)  $\lambda_{max}=492nm$ ;  $\lambda_{em}=520nm$

## Fluorophore Protein

9.0 ug/mg; 3.4 moles FITC per mole IgG

## Stabilizer

5 mg/ml BSA

## Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

## Storage

### Storage:

Store at  $-20^{\circ}C$ . Stable for one year after shipment.

### Storage Buffer:

0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 5 mg/ml BSA, 0.02% Sodium Azide, pH 7.6.

Aliquoting is unnecessary for  $-20^{\circ}C$  storage

## Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.


For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com) W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunofluorescent analysis of A431 cells, using Calnexin antibody 10427-2-AP at 1:50 dilution and SA00003-8 Fluorescein (FITC)-conjugated Donkey Anti-Rabbit IgG(H+L).