

For Research Use Only.  
Not For Use In Diagnostics.

# R-PE-conjugated F(ab')<sub>2</sub> Fragment Donkey Anti-Mouse IgG(H+L)

Catalog Number:SA00008-9



## Information

### Catalog Number:

SA00008-9

### Host:

Donkey

### Applications:

IF, FC

### Reactivity:

Mouse

### Physical State:

Liquid

### Conjugation:

R-PE

## Recommended Dilutions

1:20-1:100 for most applications

## Preservative

0.02% Sodium Azide

## Fluorophore

Phycoerythrin Amax=490nm,545nm and 566nm;Emax=580nm

## Stabilizer

5 mg/ml BSA

## Safety Notes

This product is for research use only,not for diagnostic or therapeutic use.

## Storage

### Storage:

Store at -20°C. Stable for one year after shipment.

### Storage Buffer:

0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 5 mg/ml BSA, 0.02% Sodium Azide, pH 7.6.

Aliquoting is unnecessary

## Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com) W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

