

For Research Use Only.  
Not For Use In Diagnostics.

# Cy3-conjugated Goat Anti-Rabbit IgG(H+L)

Catalog Number:SA00009-2



## Information

### Catalog Number:

SA00009-2

### Host:

Goat

### Applications:

IF, FC

### Reactivity:

Rabbit

### Physical State:

Liquid

### Conjugation:

Cy3

## Recommended Dilutions

1:20-1:100 for most applications

## Preservative

0.02% Sodium Azide

## Fluorophore

Cy3 bis-NHS ester Amax=550nm, Emax=570nm

## Stabilizer

5 mg/ml BSA

## Absorbance Ratio

A552/A280=2.4

## Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

## Storage

### Storage:

Store at -20°C. Stable for one year after shipment.

### Storage Buffer:

0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 5 mg/ml BSA, 0.02% Sodium Azide, pH 7.6.

Aliquoting is unnecessary for -20°C storage

## Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> HepG2 cells were stained with 0.2 ug Anti-Human MYH9 (11128-1-AP) and Cy3-conjugated Goat Anti-Rabbit IgG(H+L) (SA00009-2) at dilution 1:100 (red) Control Antibody in black. Cells were fixed with 4% PFA with 0.1 triton.