

For Research Use Only.  
Not For Use In Diagnostics.

# HRP-conjugated Affinipure Goat Anti-Mouse IgG, Fc $\gamma$ Subclass 1 Specific



[www.ptglab.com](http://www.ptglab.com)

Catalog Number: SA00012-1

## Information

<b>Catalog Number:</b> SA00012-1	<b>Reactivity:</b> Mouse IgG Fc $\gamma$ 1
<b>Host:</b> Goat	<b>Physical State:</b> Liquid
<b>Applications:</b> ELISA, WB, Dot blot	<b>Conjugation:</b> HRP

## Recommended Dilutions

1:2000-1:40,000 for western blotting with ECL substrates  
1:1000-1:20,000 for ELISA and Western blotting with chromogenic substrates

## Stabilizer

3 mg/ml BSA

## Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 3 mg/ml BSA, pH 7.6.  
Aliquoting is unnecessary for -20°C storage

## Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please contact:

T: 4006900926      E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)      W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data