

For Research Use Only.  
Not For Use In Diagnostics.

# HRP-conjugated Goat Anti-Mouse IgG, Fc $\gamma$ Subclass 3 Specific

Catalog Number:SA00012-5



## Information

Catalog Number:	SA00012-5	Reactivity:	Mouse IgG Fc $\gamma$ 3
Host:	Goat	Physical State:	Liquid
Applications:	ELISA, WB, Dot blot	Conjugation:	HRP

## Recommended Dilutions

1:2000-1:40,000 for western blotting with ECL substrates  
1:1000-1:20,000 for ELISA and Western blotting with chromogenic substrates

## Stabilizer

3 mg/ml BSA

## Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 3 mg/ml BSA, pH 7.6.  
Aliquoting is unnecessary for -20°C storage

## Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please contact:

T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

